



FamilyWatchToday Health Monitor

USER MANUAL

Platform: Android / Wear OS

Package: com.monitor.health

Ecosystem: FamilyWatchToday

Document Date: August 2026

Table of Contents

1.	Introduction	3
2.	Getting Started	3
	– First-Time Login	3
	– Granting Permissions	4
3.	Main Dashboard	4
4.	Health Metrics	5
	– Heart Rate	5
	– Blood Pressure	6
	– Blood Oxygen (SpO2)	6
	– Body Temperature	6
	– Blood Glucose	7
	– ECG (Electrocardiogram)	7
	– Steps	7
	– Sleep	8
	– Weight & Body Composition	8
5.	Your Profile	8
6.	Emergency & Safety Features	9
	– Automatic Fall Detection	9
	– Manual SOS	10
7.	Communication	10
	– Incoming Calls	10
	– Voice Calls	11
	– Video Calls	11
	– Messages	11
8.	Pairing Bluetooth Sensors	12
9.	Viewing Health History	12
10.	Background Monitoring	13
11.	Tips & Troubleshooting	13

1. Introduction

FamilyWatchToday Health Monitor is a wearable companion app for Android and Wear OS smartwatches. It continuously tracks your vital signs, detects emergencies such as falls, and keeps you connected with your family and care team — all from your wrist.

The app covers the following core areas:

- **Health Monitoring** — Heart rate, blood pressure, blood oxygen, temperature, glucose, ECG, steps, sleep, and weight.
- **Emergency Safety** — Automatic fall detection and one-tap SOS alerts with GPS location.
- **Communication** — Voice calls, video calls, and messages with your family or care provider.
- **Cloud Sync** — All readings are automatically saved and synced to the FamilyWatchToday platform so your care network stays informed.

Note: This app works in conjunction with the FamilyWatchToday web platform. Your account is managed by your family administrator or care provider.

2. Getting Started

2.1 First-Time Login

When you open the app for the first time, the **Login screen** will appear.

- 1 Enter your **Email address** in the first field.
- 2 Enter **Password** in the second field. Tap the eye icon on the right to your show or hide your password.
- 3 Tap **Login** . A loading indicator will appear while the app verifies your credentials.
- 4 Once confirmed, you are taken to the **Main Dashboard** .

Tip: If login fails, double-check that your email and password are typed correctly. Contact your family administrator if you do not have credentials or need a password reset.

2.2 Granting Permissions

After logging in, the app will request several device permissions. These are essential for full functionality. Tap **Allow** on each dialog.

Permission	Why It Is Needed
Body Sensors	Read heart rate, oxygen level, and other vitals from the device's built-in sensors.
Location	Send your GPS coordinates when an emergency SOS is triggered.
Bluetooth	Connect and receive readings from external BLE health sensors.
Microphone & Camera	Enable voice and video calls.
Accessibility Service	Allow the manual SOS shortcut via the home button.
Battery Optimization	Keep monitoring services running continuously in the background.

Important: If prompted to *disable battery optimization*, tap **Allow**. Without this, the app may stop monitoring when the screen is off.

3. Main Dashboard

After login, the **Main Dashboard** is the central screen of the app. It is built as a series of **swipeable pages** — each page represents one health metric.

- Swipe **left or right** to move between health metric pages.
- A brief loading screen with a spinner and the message *"Fetching data from server..."* appears while the app retrieves your latest readings.
- Once loaded, the first metric page is shown automatically.

Wear Detection

The app detects whether the watch is on your wrist. If no sensor activity is detected for 5 seconds, a **"NOT WORN"** overlay appears on screen. The overlay disappears automatically once the watch is back on your wrist and sensor contact is re-established.

Tip: Make sure the back of the watch sits snugly against your wrist for accurate sensor readings.

4. Health Metrics

The dashboard contains one page per health metric. Below is a summary of all available metrics followed by detailed descriptions of each.



Heart Rate

Real-time BPM from the wrist sensor. Tap Measure for an on-demand reading.



Blood Pressure

Systolic / Diastolic in mmHg from a paired BLE blood pressure monitor.



Blood Oxygen (SpO2)

Oxygen saturation percentage. Healthy range is typically 95–100%.



Body Temperature

Body temp in °C or °F from a paired BLE thermometer sensor.



Blood Glucose

Blood sugar level in mg/dL or mmol/L from a paired glucose meter.



ECG

Heart rhythm trace from a paired ECG sensor for cardiac monitoring.



Steps

Daily step count from the built-in accelerometer. Tap Measure to refresh.



Sleep

Duration and quality of last night's sleep from the built-in sleep tracker.



Weight & Body Composition

Weight plus estimated body fat, lean mass, water, and bone mass from a paired smart scale.

4.1 Heart Rate

The Heart Rate page shows your current **beats per minute (BPM)** with a pulsing heart animation.

Status labels:

- **Idle** — No measurement in progress.
- **Measuring** — A reading is being taken.

- **[number] BPM** — The completed reading.

Taking a Reading

- 1 Make sure the watch is properly on your wrist.
- 2 Tap the **Measure** button.
- 3 Hold still for a few seconds. The BPM value will update automatically.

Manual Entry

If you have a BPM value from another device, tap the input field, type the number, and tap **Submit** to save it.

4.2 Blood Pressure

Displays your most recent blood pressure reading as **Systolic / Diastolic** (e.g., 120 / 80) in **mmHg**, along with:

- Your pulse rate captured during the measurement.
- How long ago the reading was taken (e.g., "5 minutes ago").
- A status label: **Normal**, **High**, or **Low**.

Readings come automatically from a paired BLE blood pressure monitor (see **Section 8**).

4.3 Blood Oxygen (SpO2)

Shows your blood oxygen saturation as a **percentage** (e.g., 98%) and the pulse rate recorded during the measurement.

Health Note: A healthy SpO2 level is typically between **95% and 100%**. Contact your care provider if readings consistently fall below 95%.

Readings are captured automatically from the watch's built-in optical sensor or a paired external pulse oximeter.

4.4 Body Temperature

Displays your current body temperature in **°C** or **°F** (based on your unit preference) along with the time of the last reading.

Temperature is read from a paired external BLE thermometer and updates automatically when the sensor sends a new value.

4.5 Blood Glucose

Shows your blood glucose level in **mg/dL** or **mmol/L** along with the timestamp of the last reading.

Glucose data is received from a paired external BLE glucose meter after each measurement.

4.6 ECG (Electrocardiogram)

Displays a visual trace of your heart's electrical activity from a paired ECG sensor. The trace and timestamp are shown so you and your care provider can review your heart rhythm over time.

4.7 Steps

Shows your total **step count for the current day** from the built-in accelerometer.

Refreshing Your Step Count

- 1 Tap the **Measure** button.
- 2 The app reads the latest step count from the accelerometer.
- 3 The updated count appears on screen.

You can also type a step value manually and tap **Submit** if you need to enter a corrected count.

4.8 Sleep

Shows your sleep statistics from the previous night, including duration and quality. If your device supports sleep stage tracking, stage breakdowns (light, deep, REM) are displayed.

Tap **Measure** to sync and refresh the latest sleep data from the device sensor.

4.9 Weight & Body Composition

Shows your most recent weight in **kg** or **lbs**, along with an estimated body composition breakdown:

- **Body Fat %**
- **Lean Mass %**
- **Water %**
- **Bone Mass %**

Body composition is calculated using your weight combined with the height stored in your profile. Weight is sent automatically from a paired BLE smart scale.

Tip: To view your weight history as a chart, tap anywhere on the Weight page to open the **Weight History Graph**.

5. Your Profile

Swipe to the **Profile** page (last page on the dashboard) to see your personal information and access quick-action shortcuts.

What you see:

- Your profile photo and name.
- Personal details: age, gender, height.
- Quick action buttons (see table below).
- Current app version.

Button	Action
SOS / Emergency	Opens the emergency alert screen to send an SOS with your location.
Call	Opens the call screen.
Messages	Opens your message inbox.
Settings	Opens app and device settings.

Logging Out

To log out, open **Settings** from the Profile page and tap **Logout**.

6. Emergency & Safety Features

6.1 Automatic Fall Detection

The app continuously monitors your movement using the built-in accelerometer. If a sudden fall-like motion is detected, an emergency alert is triggered automatically.

What Happens When a Fall Is Detected

- 1

A **full-screen alert** appears with the message: *"We detected you may have fallen or requested Emergency Call."*
- 2

You have two options:

3

OK — Confirms the fall and sends an SOS alert to your family and care team.

4

Cancel — Dismisses the alert (use this if it was a false alarm).
- 5

If you tap **OK**, the following information is sent immediately:

6

Your current GPS location (latitude & longitude).

7

Device model and operating system version.

8

Current battery level.

9

Exact timestamp of the event.

False Alarms: Active exercise, dropping the device, or a sudden arm movement may occasionally trigger fall detection. Simply tap **Cancel** to dismiss — no alert will be sent.

6.2 Manual SOS

You can send an emergency alert manually at any time using either of these methods:

- **Double-tap the Home button** on your device.
- Tap the **SOS button** on the Profile page.

Both methods open the same confirmation screen as automatic fall detection. Tap **OK** to send the SOS or **Cancel** to dismiss.

Important: Tapping OK sends a real emergency alert to your family and care team. Only confirm if you genuinely need assistance.

7. Communication

7.1 Incoming Calls

When someone calls you through the FamilyWatchToday platform:

- Your screen turns on automatically — even from the lock screen.
- The caller's **name and photo** are displayed.
- A **ringtone** plays.

Button	Color	Action
Accept	Green (right side)	Answers the call.
Decline	Red (left side)	Rejects the call.

7.2 Voice Calls

Once connected on a voice call, you will see the call status and three control buttons:

Button	Action
Mute	Toggle your microphone on or off.
Speaker	Toggle speakerphone on or off.
Hang Up	End the call.

7.3 Video Calls

During a video call:

- Your **own camera feed** is shown as a preview.
- The **other person's video** fills the main screen.

All video calls use an encrypted, secure connection.

7.4 Messages

Tap **Messages** from the Profile page to open your inbox.

What You See

- Two tabs at the top: **Read** and **Unread**.

- A scrollable list of messages showing the sender name, a preview of the message, and the timestamp.

How to Use

- 1 Tap **Unread** to filter only unread messages, **Read** to see all read or ones.
- 2 Tap any message in the list to open and read the full content.

Note: Messages are sent to you by family members or care providers through the FamilyWatchToday web platform.

8. Pairing Bluetooth Sensors

The app connects to external **Bluetooth Low Energy (BLE)** health sensors for more precise and clinically validated readings. Supported sensor types:

- Blood pressure monitor
- Pulse oximeter
- Thermometer
- Glucose meter
- Smart weight scale
- ECG monitor

How to Pair a Sensor

- 1 Turn on the external sensor and **pairing mode** (refer to the sensor's own manual) put it into **mode** manual).
- 2 On the watch, go to **Settings** or the **BLE Device Selection** screen.
- 3 A list of nearby available devices will appear. Scroll to find your sensor.
- 4 Tap the sensor name.
- 5 A confirmation dialog will appear — tap **Confirm** to complete pairing.

Once paired, the sensor will automatically transmit new readings to the corresponding metric page on the dashboard.

Tip: Keep your sensors within **1–2 meters** of the watch for a reliable Bluetooth connection.

9. Viewing Health History

Weight History Graph

- 1 Navigate to the **Weight** page on the dashboard.
- 2 Tap the weight value or the graph area.
- 3 A full-screen chart opens showing your weight plotted over time with trend lines and min/max markers.

Other Metrics

Some metric pages display a mini chart or trend line below the current reading. Scroll or swipe within the page to view recent readings over time. For full historical data and detailed reports, log in to the **FamilyWatchToday web platform** from a computer or phone.

10. Background Monitoring

The app continuously monitors your health in the background — even when the watch screen is off. A small persistent notification is visible in the notification bar while monitoring is active. This is normal and expected behavior.

Service	What It Does
Heart Rate Monitor	Continuously reads BPM from the wrist sensor.
Fall Detection	Listens for fall-like motion via the accelerometer.
Step Counter	Tracks your daily step count throughout the day.
BLE Scanner	Maintains connection to paired external sensors.
Location Tracker	Periodically updates your GPS position for SOS use.
Data Sync	Uploads new readings to the FamilyWatchToday server.

Battery Optimization: Some devices stop background apps to conserve battery. If monitoring stops working after a period of inactivity, go to your device's **Settings** → **Battery** → **Apps** → **Health Monitor** and disable battery optimization for this app. The app will prompt you to do this on first launch.

11. Tips & Troubleshooting

Issue	What to Try
App shows "Fetching data from server..." for a long time	Check your Wi-Fi or mobile data. Move to an area with better signal and try again, or restart the app.
Health readings are not updating	Ensure the watch is worn snugly. For BLE sensors, confirm the sensor is powered on and within range. Tap Measure to trigger a fresh reading manually.
Fall detection keeps triggering incorrectly	Avoid sudden fast arm gestures. Tap Cancel on any false alarm — no SOS will be sent.
Cannot receive calls or messages	Make sure the app is running in the background (monitoring notification should be visible). Enable notifications for this app in device settings and check internet connectivity.
"NOT WORN" overlay will not disappear	Ensure the watch is on your wrist with the sensors against your skin. Try repositioning the watch slightly higher. The overlay should clear within a few seconds.
Bluetooth sensor not appearing in list	Make sure the sensor is powered on and in pairing mode. Move the sensor closer to the watch. Restart the sensor and try the pairing process again.
Forgot password	Contact your family administrator or care provider to reset your account credentials.
App stops monitoring after screen turns off	Disable battery optimization for the Health Monitor app in device settings (see Section 10).

For additional support, contact your **FamilyWatchToday care coordinator** or system administrator.